

Micro Quick

Ultra-Fast Particle Scanner for
Parts Cleanliness Testing



Easy to Use | Fast Results
Established Industry Standard
VDA-19.1 & ISO-16232 Compliant

by
RJL *Micro*  **Analytic**
www.particle-scanner.com

About MicroQuick

Fast and Easy to Use

- Filter analysis in 2.5 minutes plus reporting
- Front-end implemented in Microsoft Excel™

Absolute Robust

- Robust consumer technology scanner
- No warranty case within 7 years

Zero Maintenance

- Self-calibration routine to be done by user
- Automated audit report generation

VDA-19.1 & ISO-16232

The working groups of VDA-19.1 and ISO-16232 have compared the MicroQuick scanner with conventional microscopes by proficiency testing and rated the particle scanner positively.

Consequently, the scanning technology has been approved for the standard particle analysis for cleanliness testing according to VDA-19.1 (2015) and ISO-16232 (2018).

VDA | Verband der
Automobilindustrie


19

Quality Management
in the Automotive Industry


Inspection of Technical Cleanliness
– Particulate Contamination of Functionally-
Relevant Automotive Components

Three Simple Steps


① Scan



② Analysis




③ Report



Included Features


MicroEditor

- Manual refinishing (separation and stitching of particles)
- Re-classification of particles (fibres, metallic glossy)
- Length and width measure check, user rating per particle






MicroStandard

- 70 particles, size range 25-1000 µm
- 100 % detection rate guaranteed
- Automated re-calibration, reports



Metallic Gloss Detection

- Fully automatic recognition of metallic reflective particles
- Flat bed scanner equipped with unique optical filter for superior detection of metallic gloss
- RJL Micro & Analytic GmbH proprietary technology



New MicroReporter Series 5

Micro Report

Client data
Standard Components Corp.
Industrial park 3, 12345 Springfield, U.K.
Quality department, Pete Heat
Tel: +0-1234-5678-0
E-Mail: neat@standard-components.com

Work order
Reference number: 27000A/14
Quotation Number: 12345/14
Order Number: 123456789
Test objective:
Test motive:

Test object identification:
Part number:
Batch number:
Sampling:
Delivery:
Area extracted:
Number of Samples:
Surface area:
Volume:
Test directive:

Particle extraction
Analysis date:
Pre-treatment:
Technique:
Extraction device:
Nozzle type:
Nozzle Size:
Aspirator Size:
Extraction fluid:
Vol. stream:
Fluid amount:
Spray duration:
Spray angle:
Spray distance:
Feed rate:

Gravimetry
Laboratory balance:
Last calibration:
Resolution:
Measured value:

Evaluation
NOK (details see [...](#))

The results refer to:

Date: C:\... \RJLM
Report: C:\Ergebnis
Document generated

Technical Cleanliness

Testing authority
RJL Micro & Analytic GmbH
Im Entenfang 11, 76689 Karlsdorf-Neuthard
Analyseservice, Markon Graf
Tel: +49-7251-36790-0
E-Mail: zentrale@rjl-microanalytic.de

RJL Micro & Analytic
ISO 9001 DAKKS

Micro Analysis

Results

Numbers scaled	to 1 sample	Summary results	Measured	Tolerable
Numbers scaled	to 1 sample	Particle mass (mg)	1.250	1.500
Inspected lot size	1 sample	Total fibre length (mm)	36.735	
Scaling factor	1.000	Total fibre elongation (mm)	43.684	

Particle fractions

Number	Area (mm²)	Mass share estimated (mg)	Longest particle (µm)
Granules, matt	1570	17.264	1.093
Granules, metallic glossy	42	1.536	1428 x 274
Fibres, matt	89	0.935	1996 x 19
Fibres, glossy	0	0.000	-

Length-number histogram
Frequency not cumulated

Number per length class (µm)	Total	[25; 50)	[50; 100)	[100; 150)	[150; 200)	[200; 400)	[400; 600)	[600; 1000)	[1000; 1500)	[1500; 2000)	[2000; 3000)	[3000; ...)
Granules matt*	1570	527	469	218	135	172	33	13	1	1	1	0
Granules metallic glossy*	42	19	10	1	1	3	3	3	2	0	0	0
Fibres matt	89	3	6	9	8	29	16	11	5	2	0	0
Fibres glossy	0	0	0	0	0	0	0	0	0	0	0	0
Particles in selected classes*	1612	546	479	219	136	175	36	16	3	1	1	0
tolerable particles	-	-	-	-	-	-	-	20	0	0	0	0

Header

Wizard to fill-in information related to client, test object, method, parameters of extraction, filtration and instrumentation

Results

- Total particle mass and estimated per class
- Two histograms, user-defined length, width, area
- Selectable ISO-16232 or user-defined size bins
- Component cleanliness code (CCC)
- Sedimentation number for ambient monitoring
- Possibility to project numbers to specific lot size
- Length-width point diagram

Micro Gallery

Length-number
Frequency not cumulated
Number per length class (µm)
Granules matt*
Granules metallic glossy*
Fibres matt*
Fibres glossy*
Particles in selected classes*
tolerable particles
Component cleanliness
Sedimentation number

Width-length diagram

Technical Cleanliness

0) Filter membrane
Sample name: Kupplungsgehäuse

1) Granule, matt (biggest particle)
Length 2214 µm, Width 1348 µm

2) Fibre, matt
Length 1296 µm, Width 26 µm

3) Granule, matt
Length 932 µm, Width 940 µm

4) Granule, matt
Length 1469 µm, Width 917 µm

5) Granule, metallic glossy
Length 1428 µm, Width 274 µm

6) Fibre, matt
Length 1174 µm, Width 26 µm

7) Granule, metallic glossy
Length 1005 µm, Width 395 µm

8) Granule, matt
Length 997 µm, Width 637 µm

9) Granule, metallic glossy
Length 785 µm, Width 490 µm

10) Granule, metallic glossy
Length 713 µm, Width 421 µm

Gallery

- All images in colour
- Filter overview image
- Map for re-localisation of displayed particles on membrane
- Up to 100 particle micro snapshots displayed at once
- User-defined sorting by length, width, area, etc.
- Each particle can be refined by means of MicroEditor
- Selectable contrast enhancement
- Particle zoom function

Further Reports

- Mini report (single page)
- Extraction curve validation
- Cleanliness trend tool
- Scanner audit report
- ISO-16232 conformity check

MicroQuick Accessories

MicroTrap

Particle monitor for ambient cleanliness control, suitable for MicroQuick scanner



MicroStamp

Stamping tool to strip-off contamination directly from dedicated surfaces



MicroEx

Basic set of laboratory equipment for contamination extraction from parts



MicroMag

Mini microscope for manual analysis, simple navigation on membrane, magnification 50x and 200x

Multiscanner

Automated sequential particle analysis, up to three samples, suitable for:

- Filter membranes
- Particle monitor MicroTrap
- Stamping tool MicroStamp



World-Wide Customer Base



RJL Micro & Analytic GmbH
Im Entenfang 11
76689 Karlsdorf-Neuthard
Germany / Deutschland

www.rjl-microanalytic.com
sales@rjl-microanalytic.com
phone: +49-7251-36790-0
fax: +49-7251-36790-79

Sales representative